

**Accessories for Ultrasonic Homogeniser UP200St and UP200Ht**

Hielscher Ultrasonics

Accessories for Ultrasonic Homogeniser UP200St and UP200Ht

Sonotrode S26d2, titanium, approx. length 120 mm

Sonotrode S26d7, titanium, approx. length 95 mm

Sonotrode S26d14, titanium, approx. length 80 mm

Sonotrode S26d40, titanium, approx. length 80 mm

Sonotrode S26d7D, titanium, approx. length 95 mm

Sonotrode S26d26G glass (DURAN®), approx. total length 200 mm, approx. length of glass part 100 mm

Flow Cell Reactor FC7K

stainless steel flow cell reactor with cooling jacket, for sonotrodes S26d2D or S26d7D, with NBR O-rings (19 x 3,2 mm), hose connection, plastic clamps, and hose (2 m), autoclavable

Sound Protection Box SPB-L

sound protection box, acrylic glass, with divisible rod Ø12mm, vertically adjustable table, anti-slip mat, dimensions (L x W x H): 29 x 20 x 35 cm, weight: 13 kg, sound insulation: 22 dB (A) open air measurement

accessories UP200Ht/St

Code	Description	Packaging
LLG07970924	Accessories for Ultrasonic Homogeniser UP200St and UP200Ht, Type Sonotrode S26d40, Tip diam. 40 mm, For volume 100 ... 1000 ml	1 pz.
LLG07980806	Accessories for Ultrasonic Homogeniser UP200St and UP200Ht, Type Sonotrode S26d7, Tip diam. 7 mm, For volume 20 ... 500 ml	1 pz.
LLG09106625	Accessories for Ultrasonic Homogeniser UP200St and UP200Ht, Type Sonotrode S26d7D, Tip diam. 7 mm, For volume 20 ... 500 ml	1 pz.
LLG07980807	Accessories for Ultrasonic Homogeniser UP200St and UP200Ht, Type Sonotrode S26d14, Tip diam. 14 mm, For volume 50 ... 1000 ml	1 pz.
LLG07980808	Accessories for Ultrasonic Homogeniser UP200St and UP200Ht, Type Sonotrode S26d26G, Tip diam. 26 mm, For volume 50 ... 1000 ml	1 pz.
LLG07980809	Accessories for Ultrasonic Homogeniser UP200St and UP200Ht, Type Sound Protection Box SPB-L, Tip -, For volume - ml	1 pz.
LLG09106626	Accessories for Ultrasonic Homogeniser UP200St and UP200Ht, Type Flow Cell Reactor FC7K, Tip -, For volume - ml	1 pz.
LLG09106623	Accessories for Ultrasonic Homogeniser UP200St and UP200Ht, Type Sonotrode S26d2, Tip diam. 2 mm, For volume 0.15 ... 50 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

