

**Filter holders, Type 16309**

Sartorius Lab Instruments

Filter holders, Type 16309

Completely made of glass.

For particle removal from solvents. With 250 ml attachment.

For filter Ø 47 and 50 mm. With ground in vacuum bottle made of borosilicate glass, 1 l.



Code	Description	Packaging
LLG09604076	Filter holders, Type 16309, Type All-glass unit, 250 ml, Filter size Ø 47 - 50mm mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

