

**Distilling heads, DURAN tubing****Lenz****Distilling heads, DURAN® tubing**

Made of DURAN® tubing. Acc. to DIN 12594 and in addition to DIN standard (*), with one NS 14/23 socket and two NS cones (descending angle: 75°), for distillation thermometers of 55 mm built-in length.



Code	Description	Packaging
LLG06058154	Distilling heads, DURAN® tubing, Cone vertical 14/23 NS, Cone obliquely 14/23 NS, Socket 14/23 NS,	
LLG06243044	Distilling heads, DURAN® tubing, Cone vertical 24/29 NS, Cone obliquely 19/26 NS, Socket 14/23 NS, *	
LLG07613890	Distilling heads, DURAN® tubing, Cone vertical 19/26 NS, Cone obliquely 19/26 NS, Socket 14/23 NS,	
LLG09012173	Distilling heads, DURAN® tubing, Cone vertical 29/32 NS, Cone obliquely 29/32 NS, Socket 14/23 NS,	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

