

**Distilling heads, DURAN tubing**

Lenz

Distilling heads, DURAN® tubing

Made of DURAN® tubing. Acc. to DIN 12594 and in addition to DIN standard (*), with one NS 14/23 socket and two NS cones (descending angle: 75°), for distillation thermometers of 55 mm built-in length.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG06058154 | Distilling heads, DURAN® tubing, Cone vertical 14/23 NS, Cone obliquely 14/23 NS, Socket 14/23 NS, | |
| LLG06243044 | Distilling heads, DURAN® tubing, Cone vertical 24/29 NS, Cone obliquely 19/26 NS, Socket 14/23 NS, * | |
| LLG07613890 | Distilling heads, DURAN® tubing, Cone vertical 19/26 NS, Cone obliquely 19/26 NS, Socket 14/23 NS, | |
| LLG09012173 | Distilling heads, DURAN® tubing, Cone vertical 29/32 NS, Cone obliquely 29/32 NS, Socket 14/23 NS, | 1 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

