


**Distilling heads acc. to Claisen, sloping, DURAN tubing**

Lenz

**Distilling heads acc. to Claisen, sloping, DURAN® tubing**

Made of DURAN® tubing. Acc. to DIN 12594 and in addition to DIN standard (\*), with two parallel NS 14/23 sockets and two NS cones (descending angle: 75°), for distillation thermometers of 55 mm built-in length.



Code	Description	Packaging
LLG06073290	Distilling heads acc. to Claisen, sloping, DURAN® tubing, Cone vertical 24/29 NS, Cone obliquely 24/29 NS, Socket 14/23* NS	
LLG09012181	Distilling heads acc. to Claisen, sloping, DURAN® tubing, Cone vertical 14/23 NS, Cone obliquely 14/23 NS, Socket 14/23 NS	
LLG09012182	Distilling heads acc. to Claisen, sloping, DURAN® tubing, Cone vertical 29/32 NS, Cone obliquely 14/23 NS, Socket 14/23* NS	1 pz
LLG09012183	Distilling heads acc. to Claisen, sloping, DURAN® tubing, Cone vertical 29/32 NS, Cone obliquely 29/32 NS, Socket 14/23 NS	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

