



Milk sample bottles

Gerber

Milk sample bottles



| Code        | Description                                                                                 | Packaging |
|-------------|---------------------------------------------------------------------------------------------|-----------|
| LLG06266500 | Milk sample bottles, Type 50 ml, AR-glass, flat bottom, 4 ringmarks at 10 / 20 / 30 / 40 mm | 1 pz.     |
| LLG06266502 | Milk sample bottles, Type 50 ml, plastic, with snap cap                                     | 1 pz.     |
| LLG06266501 | Milk sample bottles, Type Rubber stopper for milk sample bottle 6.266 500                   | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

