



Milk sample bottles

Gerber

Milk sample bottles



Code	Description	Packaging
LLG06266500	Milk sample bottles, Type 50 ml, AR-glass, flat bottom, 4 ringmarks at 10 / 20 / 30 / 40 mm	1 pz.
LLG06266502	Milk sample bottles, Type 50 ml, plastic, with snap cap	1 pz.
LLG06266501	Milk sample bottles, Type Rubber stopper for milk sample bottle 6.266 500	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

