



Rack for Milk sample bottles, stainless steel

Gerber

Rack for Milk sample bottles, stainless steel

Meshes: 35 x 35 mm



Code	Description	Packaging
LLG06266505	Rack for Milk sample bottles, stainless steel, Type Bottlerack for 11 x 11 milk sample bottles	1 pz.
LLG06266503	Rack for Milk sample bottles, stainless steel, Type Bottlerack for 3 x 7 milk sample bottles	1 pz.
LLG06266504	Rack for Milk sample bottles, stainless steel, Type Bottlerack for 4 x 7 milk sample bottles	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

