


**Distilling link, DURAN tubing****Lenz****Distilling link, DURAN® tubing**

Made of DURAN® tubing. Distilling links without condenser, with two vertical NS-cones. With or without socket NS14/23 for distillation thermometers of 55mm built-in length.



Code	Description	Packaging
LLG09012201	Distilling link, DURAN® tubing, Effective length 150 mm, Cone 14/23 NS, Socket - NS	
LLG09012211	Distilling link, DURAN® tubing, Effective length 150 mm, Cone 14/23 NS, Socket 14/23 NS	
LLG09012213	Distilling link, DURAN® tubing, Effective length 350 mm, Cone 29/32 NS, Socket 14/23 NS	1 pz.
LLG09012203	Distilling link, DURAN® tubing, Effective length 350 mm, Cone 29/32 NS, Socket - NS	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

