




Distilling links with Liebig condenser and vacuum receiver adapter, DURAN tubing

Lenz

**Distilling links with Liebig condenser and vacuum receiver adapter,
DURAN® tubing**

Made of DURAN® tubing. DIN 12596. With Liebig condenser at 75° and vertical outlet with NS cone and vacuum connection with glass olives or thread GL14. With socket NS 14/23 for distillation thermometers of 55mm built-in length.



Code	Description	Packaging
LLG06203415	Distilling links with Liebig condenser and vacuum receiver adapter, DURAN® tubing, Effective length 160 mm, Cone 14/23 NS, Socket 14/23 NS, Connection GL	1 pz.
LLG09012261	Distilling links with Liebig condenser and vacuum receiver adapter, DURAN® tubing, Effective length 160 mm, Cone 14/23 NS, Socket 14/23 NS, Connection olive	
LLG09012263	Distilling links with Liebig condenser and vacuum receiver adapter, DURAN® tubing, Effective length 250 mm, Cone 14/23 NS, Socket 14/23 NS, Connection olive	
LLG09012264	Distilling links with Liebig condenser and vacuum receiver adapter, DURAN® tubing, Effective length 250 mm, Cone 14/23 NS, Socket 14/23 NS, Connection GL	1 pz.
LLG06203416	Distilling links with Liebig condenser and vacuum receiver adapter, DURAN® tubing, Effective length 250 mm, Cone 29/32 NS, Socket 14/23 NS, Connection GL	1 pz.
LLG09012266	Distilling links with Liebig condenser and vacuum receiver adapter, DURAN® tubing, Effective length 250 mm, Cone 29/32 NS, Socket 14/23 NS, Connection olive	1 pz.
LLG09012268	Distilling links with Liebig condenser and vacuum receiver adapter, DURAN® tubing, Effective length 400 mm, Cone 29/32 NS, Socket 14/23 NS, Connection olive	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



