

**LLG-Test and centrifuge tubes rimmed, PS or PP****LLG Labware****LLG-Test and centrifuge tubes rimmed, PS or PP**

- Made of ultra clear virgin polypropylene or virgin polystyrene
- Polystyrene tubes are transparent. Polypropylene tubes are opaque, break-resistant
- Tubes with cap are ethylene oxide sterilized
- With label

LLG-Dual-Position Caps for test- and centrifuge tubes, PE

- For culture tubes with an O.D. of 12 mm and 17 mm

2 locking positions possible:

- Closed , but gas exchange is still possible for aerobic growth conditions
- Hermetically sealed for anaerobic growth conditions, for storage or despatch of samples or for centrifugation



Code	Description	Packaging
LLG06265658	Tubes rimmed, with rim, 4 ml, Diam. 12 mm, Height 75 mm, Material PP, Max. 3000 rcf, Max. temp. 120 °C	500 pz.
LLG06265659	Tubes rimmed, with rim, 12 ml, Diam. 17 mm, Height 100 mm, Material PS, Max. 3000 rcf, Max. temp. 80 °C	250 pz.
LLG06265662	Tubes rimmed, with dual-position-cap, sterilized, 4 ml, Diam. 12 mm, Height 75 mm, Material PP, Max. 3000 rcf, Max. temp. 120 °C	500 pz.
LLG06265663	Tubes rimmed, with dual-position-cap, sterilized, 12 ml, Diam. 17 mm, Height 100 mm, Material PS, Max. 3000 rcf, Max. temp. 80 °C	500 pz.
LLG06265664	Tubes rimmed, with dual-position-cap, sterilized, 12 ml, Diam. 17 mm, Height 100 mm, Material PP, Max. 5500 rcf, Max. temp. 120 °C	500 pz.
LLG06265657	Tubes rimmed, with rim, 4 ml, Diam. 12 mm, Height 75 mm, Material PS, Max. 1400 rcf, Max. temp. 80 °C	500 pz.
LLG06265661	Tubes rimmed, with dual-position-cap, sterilized, 4 ml, Diam. 12 mm, Height 75 mm, Material PS, Max. 1400 rcf, Max. temp. 80 °C	500 pz.
LLG06265665	Tubes rimmed, dual-position-cap, ml, Diam. 12 mm, Height mm, Material PE, Max. - rcf, Max. temp. °C	1000 pz.
LLG06265666	Tubes rimmed, dual-position-cap, ml, Diam. 17 mm, Height mm, Material PE, Max. - rcf, Max. temp. °C	1000 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

