


**LLG-Transport tubes, PP, with screw cap****LLG Labware****LLG-Transport tubes, PP, with screw cap**

For convenient storage and dispatch of samples.

- Conical bottom
- Free standing
- Printed graduation
- 1 1/4 turn thread design for super fast screwing



Code	Description	Packaging
LLG06265668	Transport tubes, PP, with screw cap, 10 ml, Outer diam. 16 mm, Height with cap 80 mm	100 pz.
LLG06265667	Transport tubes, PP, with screw cap, 5 ml, Outer diam. 16 mm, Height with cap 60 mm	100 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

