




Distilling links acc. to Claisen with vacuum receiver adapter, DURAN tubing

Lenz

**Distilling links acc. to Claisen with vacuum receiver adapter, DURAN® tubing**

Made of DURAN® tubing. According to DIN 12594. With Liebig condenser (angled 75°) and vertical Claisen head, water and vacuum connection with olives. With two sockets NS 14/23 for distillation thermometers of 55mm built-in length.



Code	Description	Packaging
LLG09012271	Distilling links acc. to Claisen with vacuum receiver adapter, DURAN® tubing, Effective length 160 mm, Cone 14/23 NS, Socket 14/23 NS	1 pz.
LLG09012273	Distilling links acc. to Claisen with vacuum receiver adapter, DURAN® tubing, Effective length 250 mm, Cone 14/23 NS, Socket 14/23 NS	
LLG09012276	Distilling links acc. to Claisen with vacuum receiver adapter, DURAN® tubing, Effective length 250 mm, Cone 29/32 NS, Socket 14/23 NS	1 pz.
LLG09012278	Distilling links acc. to Claisen with vacuum receiver adapter, DURAN® tubing, Effective length 400 mm, Cone 29/32 NS, Socket 14/23 NS	



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

