


**High-performance hotplate CERAN****Gestigkeit****High-performance hotplate CERAN®**

High-performance CERAN® glass ceramic hotplate with circular heating area 145 mm diameter and square placement surface made of glass ceramic (175 x 175), stainless steel housing, power controller (adjustable 10% to 100%) and integral temperature monitor.



Code	Description	Packaging
LLG09645740	High-performance hotplate CERAN®, Type CT 10, dim. (W x D x H) 200 x 290 x 85 mm, Hotplate Ø 145 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

