

**Microfilterplates, 96-wells****AHN****Microfilterplates, 96-wells**

Glasfiber (GF) filter characteristics:

- Recommended for purification
- Prefiltering of dirty solutions
- Rapid flow rate
- High particle loading capacity

Polyethylene (PE) filter characteristics:

- Suitable for particle-removing coarse filtration
- Can be used both as a prefilter or as support filter in combination with other filters
- Hydrophobic
- Very good chemical resistance



Code	Description	Packaging
LLG06283164	Microfilterplates, 96-wells, Volume 1.0 ml, Pore size μm , Material Filter GF/FF	125 pz.
LLG06283165	Microfilterplates, 96-wells, Volume 1.0 ml, Pore size 1.6 μm , Material Filter GF/N	125 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

