



Woulff bottles, DURAN

Woulff bottles, DURAN®

- With reinforced walls, for work under vacuum
- Plastic coating to act as a splinter and implosion protection
- A glass insert is fitted with removable PP tubing connections Ø 11 mm
- Vent valve head
- Analogue pressure gauge with 2 scale ranges (1000 to 0mbar, 760mm to 0mm Hg)



Code	Description	Packaging
LLG09305341	Woulff bottles, DURAN®, Form Erlenmeyer pattern, Capacity 1000 ml, Diam. 140 mm, Neck thread 45 GL	1 pz.
LLG09305344	Woulff bottles, DURAN®, Form bottle-shaped, Capacity 10000 ml, Diam. 240 mm, Neck thread 45 GL	1 pz.
LLG09305345	Woulff bottles, DURAN®, Form bottle-shaped, Capacity 15000 ml, Diam. 255 mm, Neck thread 45 GL	1 pz.
LLG09305342	Woulff bottles, DURAN®, Form Erlenmeyer pattern, Capacity 2000 ml, Diam. 170 mm, Neck thread 45 GL	1 pz.
LLG09305346	Woulff bottles, DURAN®, Form bottle-shaped, Capacity 20000 ml, Diam. 290 mm, Neck thread 45 GL	1 pz.
LLG09305340	Woulff bottles, DURAN®, Form Erlenmeyer pattern, Capacity 500 ml, Diam. 110 mm, Neck thread 45 GL	1 pz.
LLG09305343	Woulff bottles, DURAN®, Form bottle-shaped, Capacity 5000 ml, Diam. 185 mm, Neck thread 45 GL	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

