

**Flow-through standard horn for Sonopuls Ultrasonic homogenisers****Bandelin****Flow-through standard horn for Sonopuls Ultrasonic homogenisers**

For indirect ultrasonic irradiation

- Flow-through standard horn FZ 5 G with external screwthread and titanium flat tip TT 13 FZ.
- Flow-through standard horn FZ 7 G with external screwthread and titanium flat tip TT 13 FZ.



Code	Description	Packaging
LLG09650241	Flow-through standard horn for Sonopuls Ultrasonic homogenisers, Type FZ 5 G, For HD 2070.2	1 pz.
LLG09650242	Flow-through standard horn for Sonopuls Ultrasonic homogenisers, Type FZ 7 G, For HD 2200.2	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

