

**Flow-through standard horn for Sonopuls Ultrasonic homogenisers****Bandelin****Flow-through standard horn for Sonopuls Ultrasonic homogenisers**

For indirect ultrasonic irradiation

- Flow-through standard horn FZ 5 G with external screwthread and titanium flat tip TT 13 FZ.
- Flow-through standard horn FZ 7 G with external screwthread and titanium flat tip TT 13 FZ.



| Code        | Description                                                                                 | Packaging |
|-------------|---------------------------------------------------------------------------------------------|-----------|
| LLG09650241 | Flow-through standard horn for Sonopuls Ultrasonic homogenisers, Type FZ 5 G, For HD 2070.2 | 1 pz.     |
| LLG09650242 | Flow-through standard horn for Sonopuls Ultrasonic homogenisers, Type FZ 7 G, For HD 2200.2 | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

