



Probes for Sonopuls Ultrasonic homogenisers, titanium, HD 4000 series

Probes for Sonopuls Ultrasonic homogenisers, titanium, HD 4000 series

Titanium probes.

For connection to stepped horns/booster horns.



Code	Description	Packaging
LLG06251521	Probes MS 2.0, For UW mini20, For volume 0,25 - 20 ml, For Sonpuls mini20, Diam. 2 mm, Tip microtip	1 pz.
LLG09650158	Probes TS 103, For SH 100 G SH 200 G, For volume 1 - 25 3 - 50 5 - 90 ml, For HD 4050 / 4100 / 4200, Diam. 3 mm	1 pz.
LLG09650239	Probes VS 70 T, For SH 70 G SH 213 G, For volume 10 - 200 20 - 900 ml, For HD 2070.2 / 2200.2, Diam. 13 mm, Tip probe long	1 pz.
LLG09650250	Probes TT 13, For SH 70 G SH 213 G, For volume 10 - 200 20 - 900 ml, For HD 2070.2 / 2200.2, Diam. 13 mm, Tip titanium flat tip	1 pz.
LLG09650165	Probes MS 1.5, For UW mini20, For volume 0,1 - 10 ml, For Sonpuls mini20, Diam. 1.5 mm, Tip microtip	1 pz.
LLG09650240	Probes VS 200 T, For SH 225 G, For volume 30 - 1000 500 - 2500 ml, For HD 2200.2, Diam. 25 mm, Tip probe long	1 pz.
LLG09650246	Probes KE 76, For SH 70 G SH 213 G, For volume 5 - 100 10 - 350 ml, For HD 2070.2 / 2200.2, Diam. 6 mm, Tip tapered tip	1 pz.
LLG09650243	Probes VS 190 T, For SH 219 G, For volume 25 - 900 500 - 1500 ml, For HD 2200.2, Diam. 19 mm, Tip probe long	1 pz.
LLG09650157	Probes TS 102, For SH 100 G, For volume 0,5 - 20 2 - 25 ml, For HD 4050 / 4100, Diam. 2 mm	1 pz.
LLG09650244	Probes MS 72, For SH 70 G SH 213 G, For volume 1 - 25 2 - 30 ml, For HD 2070.2 / 2200.2, Diam. 2 mm, Tip microtip	1 pz.
LLG09650249	Probes TT 25, For SH 225 G, For volume 30 - 1000 ml, For HD 2200.2, Diam. 25 mm, Tip titanium flat tip	1 pz.
LLG09650245	Probes MS 73, For SH 70 G SH 213 G, For volume 2 - 50 5 - 90 ml, For HD 2070.2 / 2200.2, Diam. 3 mm, Tip microtip	1 pz.
LLG09650159	Probes TS 104, For SH 100 G SH 200 G, For volume 3 - 50 5 - 75 5 - 100 ml, For HD 4050 / 4100 / 4200, Diam. 4.5 mm	1 pz.
LLG09650247	Probes TT 13 FZ, For FZ 5 G FZ 7 G, For volume - ml, For HD 2070.2 / 2200.2, Diam. 13 mm, Tip titanium flat tip with hole	1 pz.
LLG09650166	Probes MS 2.5, For UW mini20, For volume 0,5 - 25 ml, For Sonpuls mini20, Diam. 2.5 mm, Tip microtip	1 pz.



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

