

**LLG-Sample containers, PS, with screw cap, sterile****LLG Labware****LLG-Sample containers, PS, with screw cap, sterile**

With PE or /PP screw cap. Ideal for many applications as sample collection, transportation or storage.

- Non-cytotoxic
- EO sterilized
- 30 ml container with conical skirted bottom
- Ideal for small volume samples
- Leakproof-tested
- Unique leak tight cap ensures exceptional leak proof seal
- With plain label



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG06265646 | Sample containers, PS, with screw cap, sterile, Nominal capacity 7 ml | 700 pz |
| LLG06265647 | Sample containers, PS, with screw cap, sterile, Nominal capacity 30 ml | 50 pz |
| LLG06265648 | Sample containers, PS, with screw cap, sterile, Nominal capacity 60 ml | 60 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

