

**LLG-Sample containers, PS, with screw cap, sterile****LLG Labware****LLG-Sample containers, PS, with screw cap, sterile**

With PE or /PP screw cap. Ideal for many applications as sample collection, transportation or storage.

- Non-cytotoxic
- EO sterilized
- 30 ml container with conical skirted bottom
- Ideal for small volume samples
- Leakproof-tested
- Unique leak tight cap ensures exceptional leak proof seal
- With plain label



Code	Description	Packaging
LLG06265646	Sample containers, PS, with screw cap, sterile, Nominal capacity 7 ml	700 pz
LLG06265647	Sample containers, PS, with screw cap, sterile, Nominal capacity 30 ml	50 pz
LLG06265648	Sample containers, PS, with screw cap, sterile, Nominal capacity 60 ml	60 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

