

**LLG-Sample containers, PS, with screw cap, sterile****LLG Labware****LLG-Sample containers, PS, with screw cap, sterile**

With PE or /PP screw cap. Ideal for many applications as sample collection, transportation or storage.

- Non-cytotoxic
- EO sterilized
- 30 ml container with conical skirted bottom
- Ideal for small volume samples
- Leakproof-tested
- Unique leak tight cap ensures exceptional leak proof seal
- With plain label



Code	Description	Packaging
LLG06265646	Sample containers, PS, with screw cap, sterile, Nominal capacity 7 ml	700 pz
LLG06265647	Sample containers, PS, with screw cap, sterile, Nominal capacity 30 ml	50 pz
LLG06265648	Sample containers, PS, with screw cap, sterile, Nominal capacity 60 ml	60 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

