



## Vacuum sublimation apparatus

### Vacuum sublimation apparatus

Enables effective cleaning of chemicals and substances that turn from solids to vapours directly without liquidifying first and solidifying again on cooling. The cleaned sublimate can be fully recovered from the conical cooler without difficulty once the sublimation process is complete.



Code	Description	Packaging
LLG09305369	Vacuum sublimation apparatus, Type Macro apparatus, complete, sublimate quantity 10g to 25g	1 pz.
LLG09305367	Vacuum sublimation apparatus, Type Micro apparatus, complete, sublimate quantity 1g to 2g	1 pz.
LLG09305368	Vacuum sublimation apparatus, Type Semi-micro apparatus, complete, sublimate quantity 5g to 7g	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

