



Disposal Container SHARPSAFE, PP

Disposal Container SHARPSAFE®, PP

Sharpsafe® offers a convenient range of disposal containers and accessories that meet high requirements of hygiene and safety. The disposal of waste at its source is thus very convenient and easy.

- Container consists of 2 parts: yellow base with semi-transparent cover, safe and quick installation by simply pressing and hearing 4 clearly audible "clicks"
- The Sharpsafe® 0.6 l is a pocket container and is already pre-assembled
- The Sharpsafe® 9 l is also available with a height of 498 mm for the disposal of long items such as serological pipettes
- Equipped with a carrying handle
- Secure one-hand use possible with different mounting accessories
- Wide opening for the disposal of large objects
- Disconnection for all types of needles (2 l to 7 l), insulin pens (2 l to 7 l) and scalpel blades (2 l to 9 l)
- Opening with integrated safety flap for additional protection (2 l to 7 l)
- The closure cap can be locked in the open position if required
- 2 locking positions: provisional and final, open simply by pressing the tab
- Maximum filling level clearly visible by a black line on the lid
- Safety distance between the container opening and the clear marking line of maximum filling level
- Inscription on a label fused into the container, label will not be removed by moisture or autoclaving
- Information on the label in different languages (see table)
- The empty and open container is autoclavable up to 134 °C for 18 minutes
- Sharpsafe® is conform to ADR, TRBA 250, ISO 23907: 2013, NFX30-500 and owner of the label NF302
- Country of origin: European Union



Code	Description	Packaging
LLG09264211	Disposal Container SHARPSAFE®, PP, Type 1,0 L, 1000 ml, dim. (W x D x H) 108 x 90 x 187 mm, dim. large opening 58 x 36 mm, Labeling DE/FR/IT/ES/EN	1 pz.
LLG09264214	Disposal Container SHARPSAFE®, PP, Type 4,0 L, 4000 ml, dim. (W x D x H) 175 x 175 x 245 mm, dim. large opening 55 x 50* mm, Labeling DE/FR/IT/ES/EN	1 pz.
LLG09264213	Disposal Container SHARPSAFE®, PP, Type 3,0 L, 2600 ml, dim. (W x D x H) 197 x 120 x 205 mm, dim. large opening 55 x 50* mm, Labeling DE/FR/IT/ES/EN	1 pz.
LLG09264212	Disposal Container SHARPSAFE®, PP, Type 2,0 L, 2000 ml, dim. (W x D x H) 195 x 20 x 178 mm, dim. large opening 55 x 50* mm, Labeling DE/FR/IT/ES/EN	1 pz.
LLG09264210	Disposal Container SHARPSAFE®, PP, Type 0,6 L, 580 ml, dim. (W x D x H) 108 x 52 x 224 mm, dim. large opening 40 x 30 mm, Labeling FR/NL	1 pz.
LLG09264215	Disposal Container SHARPSAFE®, PP, Type 7,0 L, 6700 ml, dim. (W x D x H) 175 x 175 x 380 mm, dim. large opening 55 x 50* mm, Labeling DE/FR/IT/ES/EN	1 pz.
LLG09264217	Disposal Container SHARPSAFE®, PP, Type 9,0 L h, 8200 ml, dim. (W x D x H) 175 x 175 x 498 mm, dim. large opening 55 x 77 mm, Labeling FR/NL	1 pz.
LLG09264216	Disposal Container SHARPSAFE®, PP, Type 9,0 L, 9000 ml, dim. (W x D x H) 292 x 195 x 291 mm, dim. large opening 210 x 60 mm, Labeling FR/NL	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

