




O-rings

O-rings

VMQ silicone. With or without, FEP coating.



Code	Description	Packaging
LLG09142868	O-rings, Type VMQ + FEP, For flange 150 DN, Diam. 150 mm, Thickness 5 mm	1 pz
LLG09142834	O-rings, Type VMQ, For flange 60 DN, Diam. 75 mm, Thickness 4 mm	1 pz.
LLG09142846	O-rings, Type VMQ, For flange 100 DN, Diam. 110 mm, Thickness 4 mm	1 pz.
LLG09142866	O-rings, Type VMQ + FEP, For flange 60 DN, Diam. 75 mm, Thickness 4 mm	1 pz
LLG09142867	O-rings, Type VMQ + FEP, For flange 100 DN, Diam. 110 mm, Thickness 4 mm	1 pz
LLG09142857	O-rings, Type VMQ, For flange 150 DN, Diam. 150 mm, Thickness 5 mm	1 pz.
LLG09142869	O-rings, Type VMQ + FEP, For flange 200 DN, Diam. 215 mm, Thickness 5 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

