

**Magnetic stirrer Mini MR standard****IKA****Magnetic stirrer Mini MR standard**

Without heating.

- For stirring quantities up to 1L (H₂O).
- Infinitely variable speed from 0 to 2500rpm
- White top plate suitable for observing colour reactions
- Power supply with adapter for EU/CH/UK/USA

Specifications

Max. stirring volume, (H₂O): 1 l

Speed range: 0 ... 2500 rpm

Permissible ambient temperature: 5 ... 40 °C

Permissible relative moisture: 80 %

Plate material: Polyester film

Top plate dimensions (W x D): 115 x 115 mm

Overall dimensions (W x D x H): 115 x 129 x 37 mm

Weight: 0.25 kg

Supply requirements: 100 ... 240 V, 50/60 Hz

Protection to DIN EN 60529: IP42

Tested to DIN EN IEC 61010-1.



Code	Description	Packaging
LLG09720211	Magnetic stirrer Mini MR standard IKAMAG®, without heating, with EU/CH/UK/USA plug	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

