

**Magnetic stirrers, Midi MR 1 digital****IKA****Magnetic stirrers, Midi MR 1 digital**

Without heating.

- high power
- flat stable stainless steel housing
- non-locking motor
- adjustable speed
- LED digital speed display
- 0 to 56 minute timer, or continuous operation

Specifications**Motor rating input/output: 70 / 19W****Max. volume stirred (H₂O): 50L****Speed range: 0 to 1000rpm****Top plate:****Interface: 350 x 350mm, Stainless steel (AISI 304)****RS 232****Dimensions W x D x H: 360 x 430 x 110mm****Weight: 10,7kg****Protection system to DIN EN 60529: IP21****Supply requirements: 230V 50/60Hz****Tested to DIN EN IEC 61010-1.**

Code	Description	Packaging
LLG06284782	Magnetic stirrer Midi MR 1 digital without heating	1 pz.
LLG09720305	Magnetic stirrer Maxi MR 1 digital IKAMAG®, 500x500 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

