



Cell Culture Plates, treated, PS, sterile

Cell Culture Plates, treated, PS, sterile

These standard surface treated cell culture plates for the cultivation of adherent cells exhibit a reduced evaporation due to raised well rims and uniform rings on the lids.

- Available with 4-, 6-, 12-, 24-, 48-, 96- or 384-wells (lid included)
- Flat well bottom
- 96-well plate also available as round bottom plate (6.267 171)
- Single position lid reduces the risks of cross-contaminations
- Wells are labeled with alphanumeric code for easy identification
- Sterilized by gamma irradiation
- Pyrogenfree
- Material: Polystyrene



Code	Description	Packaging
LLG06267170	Cell Culture Plates, treated, PS, sterile, Working volume 0.20 ml, Culture surface 0.33 cm ² , No. of wells 96	100 pz
LLG06267171	Cell Culture Plates, treated, PS, sterile, Working volume 0.20 ml, Culture surface 0.58 cm ² , No. of wells 96	100 pz
LLG06267169	Cell Culture Plates, treated, PS, sterile, Working volume 0.29 ml, Culture surface 0.84 cm ² , No. of wells 48	100 PZ
LLG06267168	Cell Culture Plates, treated, PS, sterile, Working volume 0.57 ml, Culture surface 1.93 cm ² , No. of wells 24	
LLG06267167	Cell Culture Plates, treated, PS, sterile, Working volume 1.14 ml, Culture surface 3.85 cm ² , No. of wells 12	100 pz
LLG06267166	Cell Culture Plates, treated, PS, sterile, Working volume 2.90 ml, Culture surface 9.6 cm ² , No. of wells 6	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

