



Cell culture flasks, treated, PS, sterile

Cell culture flasks, treated, PS, sterile

These cell culture flasks with treated surfaces for the cultivation of adherent cells are available with different growth areas from 12.5 to 182 cm². Available with vented or plug-seal caps made from polyethylene. Vented caps contain a 0,22 µm hydrophobic filter to allow gas exchange and minimize risk of cross-contamination. Plug Seal caps can be used in closed systems, providing a liquid and gas sheer seal, but by simply unscrewing the cap one quarter of a turn, this cap can also be used in an open system.

Features:

- Innovative angled neck design offers good pipette and cell scraper access
- Upper triangular and wider base shape provides stability
- Stackable
- Special writing area near the neck
- Engraved graduation on both flask sides
- Sterilized by gamma irradiation
- Pyrogen-free

Material: Polystyrene



Code	Description	Packaging
LLG06267178	Cell culture flasks, treated, PS, sterile, Working volume 400 ml, Neck thread 45 mm, Culture surface 182.0 cm ² , Neck angled, Cover type standard	5 pezzi
LLG06267172	Cell culture flasks, treated, PS, sterile, Working volume 20 ml, Neck thread 45 mm, Culture surface 12.5 cm ² , Neck angled, Cover type standard	
LLG06267174	Cell culture flasks, treated, PS, sterile, Working volume 40 ml, Neck thread 45 mm, Culture surface 25.0 cm ² , Neck angled, Cover type standard	10 pz
LLG06267176	Cell culture flasks, treated, PS, sterile, Working volume 75 ml, Neck thread 45 mm, Culture surface 75.0 cm ² , Neck angled, Cover type standard	
LLG06267173	Cell culture flasks, treated, PS, sterile, Working volume 20 ml, Neck thread 45 mm, Culture surface 12.5 cm ² , Neck angled, Cover type vented	
LLG06267175	Cell culture flasks, treated, PS, sterile, Working volume 40 ml, Neck thread 45 mm, Culture surface 25.0 cm ² , Neck angled, Cover type vented	10 pz
LLG06267177	Cell culture flasks, treated, PS, sterile, Working volume 75 ml, Neck thread 45 mm, Culture surface 75.0 cm ² , Neck angled, Cover type vented	
LLG06267179	Cell culture flasks, treated, PS, sterile, Working volume 400 ml, Neck thread 45 mm, Culture surface 182.0 cm ² , Neck angled, Cover type vented	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

