

**Single channel pipette mLine, mechanical, variabel**

Sartorius Lab Instruments

**Single channel pipette mLine®, mechanical, variabel**

The mechanical pipettes of Sartorius are known for their ergonomic, user-friendly design and reliability. mLINE® pipettes represent an ergonomic choice for quality-conscious users who perform long lasting and repetitive manual pipetting. These pipettes are used in (bio)pharmaceutical laboratories, research institutes, universities, health care and industrial laboratories for various liquid handling tasks



Code	Description	Packaging
LLG06239043	Single channel pipette mLine®, mechanical, variabel, 1000 - 10000 µl	
LLG06234096	Single channel pipette mLine®, mechanical, variabel, 0.5 - 10 µl	
LLG06234095	Single channel pipette mLine®, mechanical, variabel, 0.1 - 3 µl	



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

