




Probe for temperature controllers, TEMPAT

Probe for temperature controllers, TEMPAT®

All probes with 2m connecting cable and plug.



Code	Description	Packaging
LLG09725394	Probe for temperature controllers, TEMPAT®, Type Probe Fe-CuNi, Diam. 1.5** mm, Length 250 mm	1 pz
LLG09725395	Probe for temperature controllers, TEMPAT®, Type Probe NiCr-Ni, Diam. 1.5** mm, Length 250 mm	1 pz
LLG09725320	Probe for temperature controllers, TEMPAT®, Type Probe Pt100, Diam. 1.6* mm, Length 250 mm	1 pz
LLG09725321	Probe for temperature controllers, TEMPAT®, Type Probe Pt100, Diam. 2.0* mm, Length 250 mm	1 pz
LLG09725322	Probe for temperature controllers, TEMPAT®, Type Glass-probe Pt100, Diam. 2.5 mm, Length 250 mm	1 pz
LLG09725323	Probe for temperature controllers, TEMPAT®, Type Glass-probe Pt100, Diam. 3.0 mm, Length 250 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

