

**Serial heating unit for reflow distillation****behr****Serial heating unit for reflow distillation**

Complete equipment for reflux distillation, consisting of

- Hot bench, 4 or 6 heating positions
- Standing basin
- Stand rods
- Cooling water distribution with condenser rack and cooler stand
- Reaction flask
- Behr high performance glass coolers



Code	Description	Packaging
LLG09843963	Serial heating unit for reflow distillation, Type RH 104, 4 heating places for 100ml round flasks	1 pz.
LLG09843965	Serial heating unit for reflow distillation, Type RH 254, 4 heating places for 250ml round flasks	1 pz.
LLG09843967	Serial heating unit for reflow distillation, Type RH 504, 4 heating places for 500ml round flasks	1 pz.
LLG09843964	Serial heating unit for reflow distillation, Type RH 106, 6 heating places for 100ml round flasks	1 pz.
LLG09843966	Serial heating unit for reflow distillation, Type RH 256, 6 heating places for 250ml round flasks	1 pz.
LLG09843968	Serial heating unit for reflow distillation, Type RH 506, 6 heating places for 500ml round flasks	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

