

**Vibrating shaker, VXR basic Vibrax****IKA****Vibrating shaker, VXR basic Vibrax®**

Optoelectronically controlled small shaker with a very wide speed range.

- Suitable for continuous operation
- New design and improved drive system
- Circular shaking motions
- Slow speeds are well maintained
- Attachments are interchangeable

Please order the required shaker attachment separately.

Specifications:

Shaking motion: Orbital

Orbit: 4mm

Motor rating input / output: 35 / 13.2W

Speed range:

Dimensions (W x D x H): 0 to 2200rpm

157 x 247 x 130mm

Weight (without attachment): 6.1kg

Max. loading: 2kg

Protection class DIN EN 60529: IP21

Supply requirements: 230V 50/60Hz

Tested to DIN EN IEC 61010-1.



Code	Description	Packaging
LLG06202467	Orbitalshaker Vibrax® VXR basic S 2 with UK plug	1 pz.
LLG09729900	Orbitalshaker Vibrax® VXR basic	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

