

**LLG- Beakers, low form, PTFE****LLG Labware****LLG- Beakers, low form, PTFE**

Manufactured from high quality PTFE with excellent resistance to all chemicals. Manufactured from molded, pure PTFE with a smooth internal finish. Flat base enables a flat contact to the surface for effective heat transfer. PTFE beakers have excellent resistance to chemicals and temperature stability from -200°C to +260°C.



Code	Description	Packaging
LLG06267347	Beakers, low form, PTFE, 100 ml, Diam. 55 mm, Height 64 mm	1 pz.
LLG06267348	Beakers, low form, PTFE, 250 ml, Diam. 72 mm, Height 90 mm	1 pz.
LLG06267346	Beakers, low form, PTFE, 50 ml, Diam. 45 mm, Height 55 mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

