

**LLG Syringe Filter PVDF, Polyvinylidene fluoride****LLG Labware****LLG Syringe Filter PVDF, Polyvinylidene fluoride**

We offer a line of syringe filters especially designed to provide efficient filtration of a broad variety of fluids, solvents, aqueous or inorganic solutions. Our syringe filters cover most applications in HPLC, pharmaceutical, environmental, biotechnology and food and beverage testing laboratories.

- Housing injected in Polypropylene or MABS
- Multifunctional Syringe filters equipped with Luer-Lock connections on both sides or female Luer-Lock input and male Luer-Slip output
- Suitable for all syringes with Luer connection
- Pre-sterilised or non-sterile
- Sterile products in individual hard blister
- Superior pressure stability

Hydrophilic membrane for aqueous polar and light unpolar solutions. Very low protein binding High flow rates ideal for use with biological macromolecules, water soluble oligomers and polymers. With Luer-Lock input and Luer-Slip output.



Code	Description	Packaging
LLG06285697	LLG Syringe Filter PVDF, Polyvinylidene fluoride, Type sterile, Filter diam. 13 mm, Pore size 0.2 µm, Housing Acrylic	50 pz.
LLG06285698	LLG Syringe Filter PVDF, Polyvinylidene fluoride, Type sterile, Filter diam. 13 mm, Pore size 0.45 µm, Housing Acrylic	50 pz.
LLG07970258	LLG Syringe Filter PVDF, Polyvinylidene fluoride, Type non sterile, Filter diam. 13 mm, Pore size 0.2 µm, Housing PP	500 pz
LLG09055534	LLG Syringe Filter PVDF, Polyvinylidene fluoride, Type non sterile, Filter diam. 13 mm, Pore size 0.45 µm, Housing PP	500 pz
LLG06258299	LLG Syringe Filter PVDF, Polyvinylidene fluoride, Type sterile, Filter diam. 25 mm, Pore size 0.2 µm, Housing Acrylic	50 pezzi
LLG06285702	LLG Syringe Filter PVDF, Polyvinylidene fluoride, Type sterile, Filter diam. 25 mm, Pore size 0.45 µm, Housing Acrylic	50 pezzi
LLG07970286	LLG Syringe Filter PVDF, Polyvinylidene fluoride, Type non sterile, Filter diam. 25 mm, Pore size 0.45 µm, Housing PP	500 pz
LLG07970387	LLG Syringe Filter PVDF, Polyvinylidene fluoride, Type non sterile, Filter diam. 25 mm, Pore size 0.2 µm, Housing PP	500 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

