



**Grinding Balls for PM 100 / 100 CM / 200 / 400**

**Retsch**

**Grinding Balls for PM 100 / 100 CM / 200 / 400**



| Code        | Description                                                                                    | Packaging |
|-------------|------------------------------------------------------------------------------------------------|-----------|
| LLG09738380 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Sintered aluminium oxide, Diam. 10 mm | 1 pz.     |
| LLG09738381 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Hardened steel, Diam. 10 mm           | 1 pz.     |
| LLG09738382 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel, Diam. 10 mm          | 1 pz.     |
| LLG09738383 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Agate, Diam. 10 mm                    | 1 pz.     |
| LLG09738384 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Tungsten carbide, Diam. 10 mm         | 1 pz.     |
| LLG09738400 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Silicon nitride, Diam. 10 mm          | 1 pz.     |
| LLG09738419 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, Diam. 10 mm          | 1 pz.     |
| LLG09738427 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, 0.5 kg, Diam. 2 mm   | 1 pz.     |
| LLG09738305 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Silicon nitride, Diam. 20 mm          | 1 pz.     |
| LLG09738385 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Sintered aluminium oxide, Diam. 20 mm | 1 pz.     |
| LLG09738386 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Hardened steel, Diam. 20 mm           | 1 pz.     |
| LLG09738387 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel, Diam. 20 mm          | 1 pz.     |
| LLG09738388 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Agate, Diam. 20 mm                    | 1 pz.     |
| LLG09738389 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Tungsten carbide, Diam. 20 mm         | 1 pz.     |
| LLG09738425 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel 0.5 kg, Diam. 2 mm    | 0,5 kg    |
| LLG09738426 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel 0.5 kg, Diam. 3 mm    | 0,5 kg    |
| LLG09738428 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, 0.5 kg, Diam. 3 mm   | 1 pz.     |
| LLG09738304 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Silicon nitride, Diam. 30 mm          | 1 pz.     |
| LLG09738390 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Sintered aluminium oxide, Diam. 30 mm | 1 pz.     |
| LLG09738391 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Hardened steel, Diam. 30 mm           | 1 pz.     |

| Code        | Description                                                                                    | Packaging |
|-------------|------------------------------------------------------------------------------------------------|-----------|
| LLG09738392 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel, Diam. 30 mm          | 1 pz.     |
| LLG09738393 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Agate, Diam. 30 mm                    | 1 pz.     |
| LLG09738394 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Tungsten carbide, Diam. 30 mm         | 1 pz.     |
| LLG09738417 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, Diam. 30 mm          | 1 pz.     |
| LLG09738303 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Silicon nitride, Diam. 40 mm          | 1 pz.     |
| LLG09738395 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Sintered aluminium oxide, Diam. 40 mm | 1 pz.     |
| LLG09738396 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Hardened steel, Diam. 40 mm           | 1 pz.     |
| LLG09738397 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel, Diam. 40 mm          | 1 pz.     |
| LLG09738398 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Agate, Diam. 40 mm                    | 1 pz.     |
| LLG09738399 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Tungsten carbide, Diam. 40 mm         | 1 pz.     |
| LLG09738418 | Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, Diam. 40 mm          | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

