



Grinding Balls for PM 100 / 100 CM / 200 / 400

Retsch

Grinding Balls for PM 100 / 100 CM / 200 / 400



Code	Description	Packaging
LLG09738380	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Sintered aluminium oxide, Diam. 10 mm	1 pz.
LLG09738381	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Hardened steel, Diam. 10 mm	1 pz.
LLG09738382	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel, Diam. 10 mm	1 pz.
LLG09738383	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Agate, Diam. 10 mm	1 pz.
LLG09738384	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Tungsten carbide, Diam. 10 mm	1 pz.
LLG09738400	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Silicon nitride, Diam. 10 mm	1 pz.
LLG09738419	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, Diam. 10 mm	1 pz.
LLG09738427	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, 0.5 kg, Diam. 2 mm	1 pz.
LLG09738305	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Silicon nitride, Diam. 20 mm	1 pz.
LLG09738385	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Sintered aluminium oxide, Diam. 20 mm	1 pz.
LLG09738386	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Hardened steel, Diam. 20 mm	1 pz.
LLG09738387	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel, Diam. 20 mm	1 pz.
LLG09738388	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Agate, Diam. 20 mm	1 pz.
LLG09738389	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Tungsten carbide, Diam. 20 mm	1 pz.
LLG09738425	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel 0.5 kg, Diam. 2 mm	0,5 kg
LLG09738426	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel 0.5 kg, Diam. 3 mm	0,5 kg
LLG09738428	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, 0.5 kg, Diam. 3 mm	1 pz.
LLG09738304	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Silicon nitride, Diam. 30 mm	1 pz.
LLG09738390	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Sintered aluminium oxide, Diam. 30 mm	1 pz.
LLG09738391	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Hardened steel, Diam. 30 mm	1 pz.



Code	Description	Packaging
LLG09738392	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel, Diam. 30 mm	1 pz.
LLG09738393	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Agate, Diam. 30 mm	1 pz.
LLG09738394	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Tungsten carbide, Diam. 30 mm	1 pz.
LLG09738417	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, Diam. 30 mm	1 pz.
LLG09738303	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Silicon nitride, Diam. 40 mm	1 pz.
LLG09738395	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Sintered aluminium oxide, Diam. 40 mm	1 pz.
LLG09738396	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Hardened steel, Diam. 40 mm	1 pz.
LLG09738397	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Stainless steel, Diam. 40 mm	1 pz.
LLG09738398	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Agate, Diam. 40 mm	1 pz.
LLG09738399	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Tungsten carbide, Diam. 40 mm	1 pz.
LLG09738418	Grinding Balls for PM 100 / 100 CM / 200 / 400, Material Zirconium oxide, Diam. 40 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

