

**Magnetic stirrer LLG-uniSTIRRER 1**

LLG Labware

Magnetic stirrer LLG-uniSTIRRER 1

An **ultra-flat** inductive magnetic stirrer without heating with timer. The modern magnet coil technology consists only of non-moving parts and therefore guarantees a wear-free drive.

For a better mixing of samples it is possible to change the direction of rotation every 30 seconds.

- **Timer** (1 min to 99 min, infinitely)
- 100 % wear- and maintenance-free inductive drive
- LED-display for an exact adjustment of the rotation speed
- Digital regulated rotation speed
- Alteration of the stirring direction
- Protection class IP 65
- Stainless steel surface

Scope of supply: Magnetic stirrer, power cable, adapters for UK and universal plug.

Specifications

Stirring positions: 1

Stirring quantity (H₂O): 0.8 l

Speed range: 15 ... 1500 min⁻¹

Permissible ambient temperature: 5 ... 40 °C

Magnetic bar: max. 25 mm

Display: LED

Heating plate: Ø 94 mm

Dimensions (W x D x H): 110 x 155 x 12 mm

Weight: 0.5 kg

Power supply: 100 ... 230 V, 50/60 Hz

IP code: IP 65

Warranty: 3 years



Code	Description	Packaging
LLG06263420	uniSTIRRER 1 magnetic stirrer with EU/CH, UK, AUS adapter	1 pz



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
Reagents operates with
a Certified Quality
Management System

