

**Magnetic stirrer LLG-uniSTIRRER 1****LLG Labware****Magnetic stirrer LLG-uniSTIRRER 1**

An **ultra-flat** inductive magnetic stirrer without heating with timer. The modern magnet coil technology consists only of non-moving parts and therefore guarantees a wear-free drive.

**For a better mixing of samples it is possible to change the direction of rotation every 30 seconds.**

- **Timer** (1 min to 99 min, infinitely)
- 100 % wear- and maintenance-free inductive drive
- LED-display for an exact adjustment of the rotation speed
- Digital regulated rotation speed
- Alteration of the stirring direction
- Protection class IP 65
- Stainless steel surface

**Scope of supply:** Magnetic stirrer, power cable, adapters for UK and universal plug.

**Specifications**

**Stirring positions:** 1

**Stirring quantity (H<sub>2</sub>O):** 0.8 l

**Speed range:** 15 ... 1500 min<sup>-1</sup>

**Permissible ambient temperature:** 5 ... 40 °C

**Magnetic bar:** max. 25 mm

**Display:** LED

**Heating plate:** Ø 94 mm

**Dimensions (W x D x H):** 110 x 155 x 12 mm

**Weight:** 0.5 kg

**Power supply:** 100 ... 230 V, 50/60 Hz

**IP code:** IP 65

**Warranty:** 3 years



Code	Description	Packaging
LLG06263420	uniSTIRRER 1 magnetic stirrer with EU/CH, UK, AUS adapter	1 pz



CARLO ERBA Reagents S.A.S.  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



CARLO ERBA  
Reagents operates with  
a Certified Quality  
Management System

