

**8-channel Decapper for screw cap tubes****Thermo Scientific****8-channel Decapper for screw cap tubes**

Process a full rack of Thermo Scientific Matrix or Thermo Scientific Nunc tubes in less than one minute with this ergonomic, manual capper/decapper. The Thermo Scientific 8-Channel Screw Cap Capper/Decapper is designed to cap tubes to the optimal torque to maintain sample integrity and prevent sample loss during storage. The 8-channel Capper/Decapper is a practical addition to benchtop processing of screw top storage tubes.

- Improves processing throughput of manual procedures with automated column-wise decapping and capping in just four seconds
- Ensures reproducibility by capping all tubes to the optimal torque level specified for the tube
- Allows anyone in the lab to cap or decap tubes with a light weight, easy-to-use and ergonomic design
- Stand allows easy storage of removed caps during sample processing steps
- Virtually always available for use with long battery life and quick recharge time



Code	Description	Packaging
LLG07659021	8-channel Decapper for screw cap tubes, For tubes Matrix ScrewTop (0.5 ml, 1.0 ml)	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

