



Narrow-mouth bottles with cap, PFA

Narrow-mouth bottles with cap, PFA

These unbreakable PFA bottles can be used for storage of virtually all corrosive chemicals, including hydrofluoric, nitric, and perchloric acids. PFA bottles are chemically inert, heat resistant with a non-stick, easy-to clean surface and are ideal for environmental sampling. The inert fluoropolymer cap liner ensures leak-proof sealing. Supplied with PFA caps.

- Ultra-pure, the lowest metal ion content of all fluoropolymers
- Unbreakable
- Chemically inert
- Non-stick, easy to clean
- Use for storage in pharmaceutical and semiconductor processing (intermediate stages)



Code	Description	Packaging
LLG06267438	Narrow-mouth bottles with cap, PFA, 100 ml, Neck diam. 16.7 mm, Height with cap 100,8 mm	1 pz.
LLG06267439	Narrow-mouth bottles with cap, PFA, 250 ml, Neck diam. 16.5 mm, Height with cap 133,7 mm	1 pz.
LLG06267440	Narrow-mouth bottles with cap, PFA, 500 ml, Neck diam. 17 mm, Height with cap 165,1 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

