

**Mortar grinder PULVERISETTE 2****Fritsch****Mortar grinder PULVERISETTE 2**

For wet and dry grinding of hard, medium-hard, soft, brittle and temperature-sensitive samples for analysis, quality control and materials testing. It even grinds difficult samples with a moist, fibrous or elastic structures using liquid nitrogen. And it's also excellently suited for mixing or homogenisation of organic and inorganic solids or liquids.

- Reproducible and accurate grinding results
- Precise adjustment and readout of the pestle pressure
- Short grinding and mixing times
- Loss and dust free

Max. feed size: 8 mm

Max. sample quantity: up to 190 ml

Final fineness: 10 ... 20 μm

Supply requirements: 100-120 / 200-240 V/1~, 50/60 Hz, 300 W

Dimensions (W x D x H): 310 x 460 x 410 mm

Net weight: 24 kg



Code	Description	Packaging
LLG09738200	Mortar mill "pulverisette 2", w/o.grinding tools feed volume: 30 - 150 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

