



Stirring Rods, PTFE Fluoropolymer

Stirring Rods, PTFE® Fluoropolymer

Contains fully enclosed steel insert for stiffness, yet can be bent to various shapes for use as a powered stirrer. Heat resistant to 288°C and completely inert to highly corrosive reagents. Antistick properties facilitate cleaning.

- Inert to corrosives
- Easy to clean
- Bendable into shapes



Code	Description	Packaging
LLG06267420	Stirring Rods, PTFE® Fluoropolymer, Length 204.8 mm	1 pz.
LLG06267421	Stirring Rods, PTFE® Fluoropolymer, Length 255.6 mm	1 pz
LLG06267422	Stirring Rods, PTFE® Fluoropolymer, Length 306.4 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

