


**Balls, PTFE****Saint Gobain****Balls, PTFE**

PTFE balls minimize "bumping" during boiling. They are inert and can be used up to 288°C. 100 balls per package.

- For column packing and mixing beds
- Resistant to most chemicals
- Non-contaminating



Code	Description	Packaging
LLG09012453	Balls, PTFE, Diam. 3.2 mm	100 pz.
LLG09012456	Balls, PTFE, Diam. 6.4 mm	100 pz.
LLG09012459	Balls, PTFE, Diam. 9.5 mm	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

