




Accessories for cutting mill PULVERISETTE 15

Fritsch

Accessories for cutting mill PULVERISETTE 15



Code	Description	Packaging
LLG06231269	Accessories for cutting mill PULVERISETTE 15, Type Collecting vessel, 3.5L	1 pz.
LLG09738449	Accessories for cutting mill PULVERISETTE 15, Type Collecting vessel, 60L with filter hose	1 pz.
LLG07900140	Accessories for cutting mill PULVERISETTE 15, Type Protected funnel with sample pusher	1 pz.
LLG09738424	Accessories for cutting mill PULVERISETTE 15, Type Rotor with straight cutting edges and fixed knives made of tool steel	1 pz.
LLG09738444	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, square perforation 10mm	
LLG09738441	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, square perforation 4mm	1 pz.
LLG09738442	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, square perforation 6mm	1 pz.
LLG09738443	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, square perforation 8mm	
LLG09738432	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, trapezoidal perforation 0.25mm	1 pz.
LLG09738433	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, trapezoidal perforation 0.5mm	1 pz.
LLG09738434	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, trapezoidal perforation 0.75mm	
LLG09738436	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, trapezoidal perforation 1.5mm	
LLG09738435	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, trapezoidal perforation 1mm	1 pz.
LLG09738437	Accessories for cutting mill PULVERISETTE 15, Type Sieve insert*, trapezoidal perforation 2mm	1 pz.
LLG09738429	Accessories for cutting mill PULVERISETTE 15, Type Universal support stand for the free mounting of the cutting mills	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

