

**Coil condensers, DURAN tubing****Lenz****Coil condensers, DURAN® tubing**

Made of DURAN® tubing. Acc. to DIN 12593 and in addition to DIN standard.  
With jacketed cooling coil. Water connection with glass olives or GL 14 thread  
with two threaded plastic olives.



Code	Description	Packaging
LLG06304794	Coil condensers, DURAN® tubing, Effective length 160 mm, Cone 14/23 NS, Socket 14/23 NS, Connection olives, DIN, Cooling surface 240 cm²	
LLG06227832	Coil condensers, DURAN® tubing, Effective length 250 mm, Cone 29/32 NS, Socket 29/32 NS, Connection olives, DIN 12593, Cooling surface 400 cm²	1 pz.
LLG06240342	Coil condensers, DURAN® tubing, Effective length 400 mm, Cone 29/32 NS, Socket 29/32 NS, Connection olives, DIN 12593, Cooling surface 650 cm²	1 pz.
LLG07635816	Coil condensers, DURAN® tubing, Effective length 160 mm, Cone 29/32 NS, Socket 29/32 NS, Connection olives, DIN 12593, Cooling surface 240 cm²	1 pz.
LLG09012542	Coil condensers, DURAN® tubing, Effective length 160 mm, Cone 29/32 NS, Socket 29/32 NS, Connection GL14, DIN 12593, Cooling surface 240 cm²	1 pz.
LLG09012543	Coil condensers, DURAN® tubing, Effective length 250 mm, Cone 29/32 NS, Socket 29/32 NS, Connection GL14, DIN 12593, Cooling surface 400 cm²	1 pz.
LLG09012544	Coil condensers, DURAN® tubing, Effective length 400 mm, Cone 29/32 NS, Socket 29/32 NS, Connection GL14, DIN 12593, Cooling surface 650 cm²	1 pz.
LLG07970241	Coil condensers, DURAN® tubing, Effective length 500 mm, Cone 29/32 NS, Socket 29/32 NS, Connection olives, DIN, Cooling surface 840 cm²	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

