

**Dipper/Sieve Chemware, PFA****Saint Gobain****Dipper/Sieve Chemware®, PFA**

Resistant to most chemicals, the articles are ideal for ultra-pure applications and use with aggressive media. Handle design allows dipper to hook over edge of vessel, eliminating dripping of corrosive fluids onto surrounding surface. Use with containers that have at least an 80 mm opening. The sieve mesh size is available in 1.0 mm.

- Non-contaminating and resistant to most chemicals
- Superior chemical resistance
- Lowest concentration of metal ions of any fluoropolymer
- Ideal for coarse sorting or draining of aggressive media
- Handle hooks over edge of vessel for draining



Code	Description	Packaging
LLG06267475	Dipper/Sieve Chemware®, PFA, Diam. 76 mm, Length 180 mm, Dipper	1 pz.
LLG06267477	Dipper Chemware®, mw.1,0mm PFA, dia. 76	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

