

**pH electrode Eutech, not refillable****Thermo Scientific****pH electrode Eutechä, not refillable**

General purpose plastic-body pH electrode. Used with all pH meters with BNC input connector.

Specifications

Temperature range: 0 ... 80 °C

pH range: 0 ... 13

Temperature sensor: no

Length / Diam.: 90 x 12 mm

Electrolyte: Reference Electrolyte

Diaphragm: Porous HDPE pin



Code	Description	Packaging
LLG06267476	Dipper Chemware®, mw.0,5mm PFA, dia. 76	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

