

**Cleanroom Wipes Sontara MicroPure, polyester/cellulose****Sontara****Cleanroom Wipes Sontara® MicroPure, polyester/cellulose**

Sontara® MicroPure nonwoven styles (woodpulp/polyester) are engineered for cleanroom use and combine absorbency of a natural fibre with the cleanliness and strength of a synthetic fibre in the same fabric. No binders or chemical treatments are used to manufacture them. The high cleanliness and purity are ensured by low extractable level and ion content as well as extremely low particle generation by the material itself. The fabric styles offer high particle and liquid absorbency due to their cellulose content.



Code	Description	Packaging
LLG09390185	Cleanroom Wipes Sontara® MicroPure, polyester/cellulose, Type Sontara® MicroPure 100, dim. 305 x 305 mm, Cleanroom Class ISO 5 + 6	150 pz.
LLG09390186	Cleanroom Wipes Sontara® MicroPure, polyester/cellulose, Type Sontara® MicroPure 100, dim. 224 x 224 mm, Cleanroom Class ISO 5 + 6	300 pz.
LLG09390187	Cleanroom Wipes Sontara® MicroPure, polyester/cellulose, Type Sontara® MicroPure AP, dim. 305 x 305 mm, Cleanroom Class ISO 6 - 9	150 pz.
LLG09390188	Cleanroom Wipes Sontara® MicroPure, polyester/cellulose, Type Sontara® MicroPure AP, dim. 224 x 224 mm, Cleanroom Class ISO 6 - 9	300 pz.
LLG09390189	Cleanroom Wipes Sontara® MicroPure, polyester/cellulose, Type Sontara® MicroPure SV, dim. 305 x 305 mm, Cleanroom Class ISO 7 - 9	250 pz.
LLG09390190	Cleanroom Wipes Sontara® MicroPure, polyester/cellulose, Type Sontara® MicroPure SV, dim. 224 x 224 mm, Cleanroom Class ISO 7 - 9	500 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

