

**LLG-Cleaning brushes****LLG Labware****LLG-Cleaning brushes**

Excellent tools for washing and cleaning of all type of laboratory bottles, flasks, tubes, condensers, burettes, etc. Brush material is manufactured from highest quality of fibre/perlon for long life cycle. Brush stem is manufactured from corrosion proof twisted steel wire. Brush handle is manufactured from polypropylene and offers extreme easy handling and use of the brush. Unique handle design enables hanging of the brush to a convenient place for easy reach. Spherical sponge top part enables access to every corner of the labware and guarantees perfect cleaning results.



Code	Description	Packaging
LLG06267298	Cleaning brushes, For narrow test tubes, Cleaning brush 12x80mm, Length 235 mm, Figure f	10 pz.
LLG06267294	Cleaning brushes, For large size bottles, Cleaning brush 70x130mm, Length 460 mm, Figure b	10 pz
LLG06267296	Cleaning brushes, For small size bottles, Cleaning brush 40x110mm, Length 380 mm, Figure d	10 pezzi
LLG06267293	Cleaning brushes, For burettes and tubings, Cleaning brush 20x125mm, Length 845 mm, Figure a	10 pz.
LLG06267297	Cleaning brushes, For test tubes, Cleaning brush 25x100mm, Length 295 mm, Figure e	10 pezzi
LLG06267295	Cleaning brushes, For medium size bottles, Cleaning brush 50x110mm, Length 460 mm, Figure c	10 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

